

## Creating the Platform

Awhile back, physical therapists made us aware that they needed a way to connect bands and various other therapeutic innovations to their patients via our harness. We designed a system to attach bands and other assistive devices, so our harness evolves as a **virtual platform** for physical therapy.

You're looking at the beginning of our **PT Adapter Kit**. We are offering these unique clips, rings, cuffs and connectors as a way to utilize your innovations through our harness - to aid in rehabilitation.

Our unique fasteners connect into the existing clips of the harness, thereby providing additional clips and rings to enable **PT's** to connect their choice of therapeutic bands or assistive devices onto the harness. Additional handles and future innovations like custom ice packs, water treadmill connectors and dog cart connector straps are coming.

All of it starts with this, a simple way to connect virtually anything to the harness. We're hoping you will take these innovations and run with them. Just let us know what you need, **so that we can help you help them!**

**Resistance Bands not included.**

**Also watch our fitting demonstration at [www.helpemup.com/videos](http://www.helpemup.com/videos)**

**Blue Dog Designs 1550 South Acoma Street, Denver CO 80223 720-237-6852 [www.helpemup.com](http://www.helpemup.com)**



Gracie, 8

## Introducing The **PT Adapter Kit** Therapeutic Exercise Connectors



Help them by helping them **up**.

Charlotte, 9

# THE HELP **EM** UP HARNESS

The lifting & mobility harness recommended by most **Vets & Therapists.**

## Multi-Clips with "O" Rings

Sizes:



**Multi-Clips with "O" Rings** These clips allow you to connect any available resistance strap by looping them through the rings or connecting it to them. These Multi rings clip into all harnesses clips.

## Connector Cuffs with "O" Rings

Sizes:



Good for Large to Extra Large sized dogs  
65 - 200 lbs

Good for Small to Medium sized dogs  
10 - 65 lbs



**Connector Cuffs with "O" Rings** have a soft textured interior pad that inhibits them from sliding. They also utilize a Velcro tightening connection that lets you easily adjust to body locations.

## Connector Clip with "O" Rings

Size:



This **Connector Clip with "O" Rings** allows you to quickly attach devices to any part of our harness, to "D" rings, "O" rings or handles. Allowing you to change angles to obtain optimal resistance.



**Right:** A Multi Clip with O Rings is connected to the belly band - and a Biko Band runs forward to a Medium Cuff.

**Below:** Once the Multi Clip is connected to the chest, reconnect the strap to the clip under the rings.

**Below right:** Future water tank straps will allow you to support, center and release dogs much more easily.



**Right:** A Connector Cuff with "O" Rings is wrapped around the hind foot just below the hock.

**Below:** Shown in detail, Connector Cuffs with "O" Rings let you adjust tension very easily. The dual rings allow you to tighten or relax the band to add or reduce tension.



**Right:** Angles, it's all about angles! This clip lets you find the perfect angle to connect therapy bands for each pup. It's also works for attaching other assistive devices.

**Below:** Most of you will want to connect two at a time for balanced resistance. These clips make it easier to apply to a moving dog.



**Above:** To apply one, first disconnect the strap on the chest or hip lift, then plug the Multi-Rings back into the connector and then reconnect the strap back into the receptor clip.

**Left:** Shown in detail, Multi-clips with "O" Rings at the chest clip with resistive band linked to Connector Cuff on rear leg.



**Multi-Clips with "O" Rings** allow you to connect virtually all resistance bands - especially when used with our Connector Cuffs. We also know you're going to find many other ways to apply them.

**Multi-Clips with "O" Rings** will have future uses that will include looped handles at front and back or across the hips. Also coming will be connecting straps for carts, water therapy tanks and ice packs.



**Connector Cuffs** let you connect to the front and back legs or side to side to add resistance. These cuffs can also help in the future with additional innovations.



**Resistance Bands not included.** Bands are sold separately and are not marketed by Blue Dog Designs.