



Size	Pad Width	Adjusting	Weight*	Gender	
Mini	4 cm	71 - 175 cm	Under 7 kg	FEMALE	MALE
Small Male	9 cm	94 - 152 cm	7 - 20 kg		
Small Female	10 cm	94 - 152 cm	7 - 20 kg	Carectellan	Company Inc.
Medium (Large Male)	14 cm	96.5 - 152 cm	20 - 90 kg		
Large Female	19 cm	96.5 - 152 cm	25 - 90 kg		400
Tall Male (For Lean Dogs)	12.5 cm	66 - 94 cm	32 - 90 kg		
Tall Female (For Lean Dogs)	19 cm	66 - 94 cm	32 - 90 kg		

### \*Size should be determined by Pad Width and Adjusting Measurements, not the weight of your dog

# How to Measure

# **Support Pad Width**

Measure from the front of the hind legs (in front of the penis for male dogs) towards your dogs chest to determine the support pad width. For dogs with barrel chests (i.e. Boxers, Pit Bulls), measure just the belly area to the base of the rib cage. Downsizing is common for these dogs.



# chest harness not included

# **Adjusting Dimensions**

Measure starting where the handle of the GingerLead will rest in your hand (while you are standing upright), continue under your dog's belly, up the other side, then back to your hand completing the loop. If it measures 94cm or less, the Tall size should be considered.

# Still not sure? ...try the towel test

## **Female Dogs**

Fold a bath towel to 19cm wide to simulate a Large Female GingerLead, then sling it under her belly. If that's too big, then re-fold the towel to 10cm wide to simulate a Small Female GingerLead. If she is in between sizes, re-fold the towel to 14cm to simulate a Medium (Large Male) GingerLead. If the small is too big for her, try our Mini GingerLead.

## Male Dogs

Fold a bath towel to 14cm wide to simulate a Medium (Large Male) GingerLead, then sling it under your dog's belly in front of his penis. If that's too big, then re-fold the towel to 9cm wide to simulate a Small Male GingerLead. If that's too big for him, try our Mini GingerLead. Male Dachshunds typically use the Mini GingerLead with the sling behind their penis.

Check with your veterinarian if you have any questions about whether a GingerLead will fit your dog's needs given their specific medical condition(s).